Integrity® Multi•Cure Temporary Crown & Bridge Material

DIRECTIONS FOR USE – ENGLISH

For Dental Use Only
USA: Rx Only

1. PRODUCT DESCRIPTION

Integrity® Multi•Cure Temporary Crown and Bridge Material is an automixed two-component material based on multi-functional methacrylic esters, intended for use as both a short- and long-term temporary crown and bridge material. Integrity® Multi•Cure Temporary Crown and Bridge Material contains both autocure initiators and photoinitiators, providing reliable self curing of the material, with the operator's option of utilizing a light curing unit to finalize the cure, reducing treatment time.

1.1 Delivery forms

Integrity® Multi•Cure Temporary Crown and Bridge Material is available in:

- 76g Cartridges for use with 10:1 automixing gun
- 7 shades – BW, B1, A1, A2, A3, A3.5, C2

1.2 Composition

Glass filler; Methacrylate Monomers; Catalyst; Photoinitiator; Stabilizers

1.3 Indications

- Temporary crowns, bridges, inlays, onlays, veneers
- Long-term temporaries

1.4 Contraindications

Integrity® Multi•Cure Temporary Crown and Bridge Material is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins.

1.5 Compatible temporary cements

Integrity® Multi•Cure Temporary Crown and Bridge Material temporary restorations can be cemented with most commonly used temporary luting agents. See manufacturer's complete Directions for Use for compatibility and specific technique. Most eugenol-containing temporary cements are compatible with Integrity® Multi•Cure Temporary Crown and Bridge Material. However, they can have a deleterious influence on resin-containing luting systems chosen for final restoration cementation.

1.6 Compatible restorative materials

Integrity® Multi•Cure Temporary Crown and Bridge Material is chemically compatible with conventional methacrylate-based composite restorative materials for marginal repair and aesthetic contouring and/or custom shading of temporary crowns. See manufacturer's complete Directions for Use for compatibility and specific technique.

1.7 Technical specifications and general material timing

For intraoral fabrication/extraoral curing, finishing, polishing:

- Immediately before each use. Bleed cartridge prior to inserting new mixing tip
- 0-15 seconds Material dispensed into preliminary impression
- 15-45 seconds Material inserted into mouth
- 1'30"-2'15" Removal from mouth

Restoration Curing:

- LIGHT CURE- spectral output containing 450nm, minimum light output at least 400mW/cm²
- 20 seconds per unit with hand-held curing light **OR** 60 seconds in extraoral VLC curing oven

CHEMICAL (SELF) CURE

- 5 minutes

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other sections of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Integrity® Multi•Cure Temporary Crown and Bridge Material contains polymerizable methacrylate monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
 - **Avoid Eye Contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - **Avoid Skin Contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and alcohol and wash well thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - **Avoid Contact with Oral Soft Tissues/Mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
2. The polymerization reaction of Integrity® Multi•Cure Temporary Crown and Bridge Material is exothermic. Do not apply directly to pulp tissue. Do not light cure material directly on or against tooth structure to prevent overheating of tooth preparation (see Step-By-Step Instructions).

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.
2. Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
3. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
4. Material should extrude easily: **DO NOT USE EXCESSIVE FORCE**. Excessive pressure may result in unanticipated extrusion of the material or cause cartridge rupture.

5. Interactions:

- Integrity® Multi•Cure Temporary Crown and Bridge Material will adhesively bond to most dental adhesives. Use of a dentin/enamel adhesive on the preparation prior to restoration fabrication may make restoration removal and cleanup difficult.
- Integrity® Multi•Cure Temporary Crown and Bridge Material may adhesively bond to freshly placed methacrylate and composite core buildup materials. Buildups should be fully cured, prepared, and rinsed thoroughly with forceful air/water spray before application of Integrity® Multi•Cure Temporary Crown and Bridge Material. If surfaces were not prepared or instrumented, wipe thoroughly with alcohol soaked gauze, then rinse with air/water spray.
- Residual Integrity® Multi•Cure Temporary Crown and Bridge Material may interfere with the setting of some Vinyl Polysiloxane impression or bite registration materials. After initial restoration fabrication, prior to impression/registration, tooth preparations should be wiped thoroughly with alcohol soaked gauze, then rinsed with air/water spray.
- Eugenol- and hydrogen peroxide-containing materials should not be used in conjunction with this product since they may interfere with hardening and cause softening of the polymeric components of the material.
- Integrity® Multi•Cure Temporary Crown and Bridge Material is a light-cured material. Therefore, it should be protected from ambient light.

2.3 Adverse reactions

- Product may irritate the eyes and skin. **Eye contact:** irritation and possible corneal damage. **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin. **Mucous membranes:** inflammation (see Warnings).
- Product may cause pulpal effects (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Keep out of direct sunlight and store in a well ventilated place at temperatures between 15-25°C/59-77°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS FOR USE

Cartridge Dispensing Gun Instructions

1. Insert cartridge into dispensing gun
2. Bleed cartridge before installing new mixing tip
3. Dispensing
 - Squeeze dispenser trigger with moderate and even pressure to begin mixing/dispensing material. Dispense a small amount through the mixing tip onto a mixing pad and discard. Proceed immediately to application.
4. Disinfect cartridge dispensing gun, cartridge and attached used mixing tip for storage.

Temporary Restoration Instructions

1. Preliminary impression
 - 1.1 Make an alginate or silicone impression of the quadrant including tooth/teeth to be prepared following manufacturer's Directions for Use.
 - **Technique Tip:** Alternatively, a vacuum-formed matrix may be constructed on a cast (with diagnostic wax-up, contour modification as indicated).
 - 1.2 Trim the interdental gingival areas out of the preliminary impression in the area(s) to be restored to provide bulk for the later finishing of the Integrity® Multi•Cure Temporary Crown and Bridge restoration.
2. Prepare tooth/teeth
 - 2.1 Ensure adequate reduction for material thickness, especially on functional occlusal surfaces.
 - 2.2 Trial-fit impression (or vacuum-formed matrix) made in Step 1 to ensure adequate orientation. Rinse impression and lightly air dry.
 - 2.3 Block out preparation undercuts with wax or suitable material, clean tooth/teeth preparation with alcohol soaked gauze and/or air water spray and lightly dry. Do not desiccate tooth structure.
 - **Technique Tip:** A thin layer of water-soluble lubricating medium may be applied to the tooth preparation to facilitate removal of the curing Integrity® Multi•Cure Temporary Crown and Bridge material.
3. Mixing and application
 - **NOTE:** Allow Integrity® Multi•Cure Temporary Crown and Bridge Material to reach room temperature prior to use.
 - 3.1 Follow instructions for cartridge loading, bleeding and mixing tip installation. Use only the supplied mixing tips.
 - 3.2 Dispense a small amount through the mixing tip onto a mixing pad and discard. Without delay, begin dispensing Integrity® Multi•Cure Temporary Crown and Bridge material directly into the preliminary impression into the area of the preparation by injecting it into the deepest portion of the preliminary impression and then into the gingival areas, keeping the mixing tip immersed in the material to prevent bubbles. Reseat the preliminary impression into the patient's mouth in the proper relative position within 45 seconds of the start of mix.
 - 3.3 Leave impression in place for approximately 1'30".
4. Removal from mouth
 - 4.1 Remove the impression containing the Integrity® Multi•Cure Temporary Crown and Bridge restoration during the firm elastic stage. Monitor the curing process carefully because removal of the temporary restoration is only possible during the elastic phase. **NOTE:** Mouth temperature significantly influences set time. If the material remains on the preparation when the impression is removed, gently tease the restoration off with gauze. Use finger pressure only, and dislodge along the path of insertion. Do not rock or torque restoration off the preparation. Avoid the use of hemostats or other pliers that may distort the setting Integrity® Multi•Cure Temporary Crown and Bridge restoration.
 - 4.2 At this elastic phase, gross excess may be quickly trimmed with sharp scissors, being careful not to distort the restoration or margins.
5. Curing
 - 5.1 Do not replace restoration onto prepared tooth/teeth for final curing. Allow restoration to remain undisturbed, within or out of the matrix, as removed from the mouth, until complete set, approximately 5 minutes from start of mix. Alternatively, light cure the restoration with a curing light for 0'20" per unit or place impression containing restoration in an extraoral light curing oven for 0'60".
6. Finishing and polishing
 - 6.1 Remove the excess material and proximal undercuts with rotary instruments.
 - 6.2 Restoration may be polished with flour of pumice and acrylic polish or by use of Enhance® Finishing System and PoGo® Polishing System or Prisma® Gloss® Composite Polishing Paste System. (See complete Directions for Use supplied with polishing product chosen)
7. Luting the temporary restoration
 - 7.1 Check the fit and esthetics of the Integrity® Multi•Cure Temporary Crown and Bridge provisional restoration. If additions are necessary, see Repairing/Modifying the Temporary Restoration. The internal surface of the restoration should be clean and dry prior to cementation.
 - 7.2 Follow temporary luting agent manufacturer's complete Directions for Use for specific technique.

Repairing/Modifying the Temporary Restoration

1. Freshly made restorations

For adjustment of margin discrepancies, roughen the area, extending slightly beyond the margin onto restoration surface. Apply a fresh mix of Integrity Multi•Cure Temporary Crown and Bridge material or a compatible flowable composite and light cure as outlined in the complete Directions for Use. For surface contour and/or esthetic adjustment, cut back area to be modified with a coarse diamond. Apply and cure desired shade(s) of characterization material as outlined in the complete Directions for Use. Adhesive/Bonding agent application is not necessary for modifications to Integrity Multi•Cure temporary crown and bridge restorations. Finish, polish and cement as outlined above.

2. Restorations in service

Integrity Multi•Cure temporary crown and bridge restorations are characterized by excellent mechanical properties. However, should breakage occur, the following procedures are recommended.

- Roughen the surfaces of the breakpoint with a coarse diamond bur and place undercuts in the adjacent areas. Prepared parts may be placed intraorally to assure fit and proper alignment. If fractured sections do not approximate, consideration should be given to remaking, rather than repairing the restoration.
- Use freshly mixed Integrity Multi•Cure Temporary Crown and Bridge Material to connect the parts. Light cure the repair area for 0'20" or hold in place for 1 minute and allow to fully cure for a minimum of 5 minutes. Remove and finish and polish as outlined above.

Technique Tips to avoid breakage:

- Be sure adequate preparation reduction/occlusal clearance exists
- Adjust occlusion removing undesirable contacts
- Allow adequate bulk for connectors in multiple unit provisional restorations
- Avoid squeezing the Integrity Multi•Cure Temporary Crown and Bridge restoration with pliers when removing restoration
- Integrity Multi•Cure Temporary Crown and Bridge restorations are durable under usual provisional conditions. Many factors contribute to the longevity, thus no definitive lifespan can be stated. However, if final restorations are not planned to follow in timely succession, consideration must be given to subsequent removal and restoration remake.

4. HYGIENE

Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

Cleaning the dispenser

The cartridge dispensing guns may be cleaned by scrubbing with hot water and soap or detergent. Excess material may be cleaned with alcohol-moistened gauze.

Disinfection of the dispenser gun and cartridge

Disinfect dispenser guns, cartridges and used tips left in place on cartridges for storage with a hospital-level tuberculocidal disinfectant according to national and local regulations. Use impregnated wipes not sprays. **NOTE:** Vigorous wiping can destroy the label. Wipe cartridge gently.

Some phenolic based agents and iodophor based products may cause surface staining. Agents containing organic solvents, such as alcohol may tend to dissolve the plastic dispenser and plunger. The disinfectant manufacturer's directions should be followed properly for optimum results.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard uses: "YYYY/MM."
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date

 Manufactured by
Dentsply Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
Made in USA
1-302-422-4511
www.dentsplysirona.com

Dentsply DeTrey GmbH
De-Trey-Strasse 1
78467 Konstanz
Germany
49-7531-583-0
www.dentsplysirona.com

Distributed by:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada

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